

alinak@bgu.ac.il linakarabchevsky.com


Glowing Microfluidics

A. Karabchevsky^{1*}, A. Mosayyebi², A. V. Kavokin^{3,4}

- ¹ Electrooptical Engineering, Ben-Gurion University, Beer-Sheva, Israel
- ² Optoelectronics Research Centre, University of Southampton, UK
- ³ School of Physics and Astronomy, University of Southampton, UK
- ⁴ CNR-SPIN, Viale del Politecnico 1, I-00133 Rome, Italy

Introduction


Chemiluminescence is a fascinating optical finds its use in application areas: from forensic science to industrial bio-chemistry. Luminol is chemical that exhibits chemiluminescence with emitted blue glow light.



The chemiluminescence reaction of luminol C₈H₇N₃O₂ (I) in relatively nonacidic solvents resulting in aminophthalate ion (III) which is a light emitting species


Materials and Methods

We have designed and fabricated a reusable microflow^{4,5} device with a serpentine channel of 600 mm width, 200 mm depth and 600 mm



Results

Here we observe an enhancement of the chemiluminescence intensity of a luminol flow metallic the presence of nanospheres on a microfluidics Experimental results below captured by the CCD camera (a) shows glowing serpentine channel compared to the Reference. Luminol has been injected with the flow rate of 0.35 mL/sec.



Conclusions

To conclude, we proposed two possible for chemiluminescence mechanisms enhancement during flow injection: 1) The rate of emission by a chemophore emitter can be enhanced by an antenna effect. The Purcell effect may be responsible for the amplification of the radiative decay rate by the enhanced density of optical states accessible for decay of the molecular excitation. 2) Observed enhancement and modification of the spectral shape could be also due to the far-field scattering.

References

1. Karabchevsky A, Mosayyebi A, Kavokin AK. Tuning chmiluminiscence of a luminol flow by plasmonic nanoparticles. Light Science and Applications; Under Review.

3

01 405 435 465 495 525 555 585 615 645 Wavelength (nm)

30 nm

gold and silver $\frac{3}{2}$ 0.2

nanoparticles

acceptors.

- 2. Agio M, Cano DM. Nano-optics: the Purcell factor of nanoresonators. Nat Photon 2013; 7: 674-675.
- 3. Greffet JJ. Nanoantennas for light emission. Science 2005; 308: 1561-1563.
- 4. Bhattacharyya A, Klapperich CM. Design and testing of a disposable microfluidic chemiluminescent immunoassay for disease biomarkers in human serum samples. Biomed Microdevices 2007: 9: 245-251.
- 5. Sackmann EK, Fulton AL, Beebe DJ. The present and future role of microfluidics in biomedical research. Nature 2014; 507: 181-189.

